

Supplementary Table S1. Coefficients of variability (CV) for calprotectin (CP) and myeloperoxidase (MPO) ELISAs. Inter-assay CV were calculated from each sample (high, medium and low) assayed in duplicate over 5 days. Intra-assay CV were calculated from repeated analysis of each sample (high, medium and low) 10 times on the same plate.

<i>Sample</i>	<i>Inter-assay CV</i>		<i>Intra-assay CV</i>	
	<i>CP</i>	<i>MPO</i>	<i>CP</i>	<i>MPO</i>
<i>High</i>	54.8%	19.9%	10.4%	0.8%
<i>Medium</i>	62.5%	NA	31.4%	34.8%
<i>Low</i>	58.7%	47.0%	14.2%	34.5%